MEDITAPE™ UC-12S

Identification of the IVD reagent MEDITAPE™ UC-12S

Intended use

For in vitro diagnostic use only
MEDITAPE UC-12S is a urinalysis test strip with reagent pads for the determination of diagnostic parameters in human urine.
MEDITAPE UC-12S test strips are for visual measurement and for use with the following Sysmex Urine Chemistry Analyzer only:
UC-1000

Principles of the examination method¹⁾

[Urobilinogen] Azo coupling method: The reaction of urobilinogen in urine with 3,4-Methylene dioxybenzenediazonium tetrafluoroborate develops a color change from light pink to red.
[Blood] Peroxidase-like reaction of the hemoglobin: Peroxidase-like hemoglobin releases oxygen, oxidizing tetramethylbenzidine causing the color change from white to blue.
[Protein] Protein error of pH indicator: A reaction between tetrabromophenol blue and protein induces a color change from yellow to blue-green.
[Glucose] Enzyme method (GOD, POD): Glucose in the urine is determined by the glucose oxidase-peroxidase reaction. The action of glucose oxidase on glucose produces gluconic acid and hydrogen peroxide. The peroxidase then catalyzed the reaction of hydrogen peroxide and chromogen to form a blue dye. The basic color on the test pad is yellow and the color will change from yellow to blue-green according to the glucose concentration.
[Ketones] Nitroprusside method: Acetone and acetoacetic acid react with sodium nitroprusside in an alkaline medium. A color change from buff pink to lavender occurs.
[Bilirubin] Azo coupling method: 2,4-dichlorobenzene diazonium tetrafluoroborate reacts with bilirubin in the presence of acid medium. A color change from light tan to light pink occurs.
[Nitrite] Griess method: Nitrite forms a diazonium compound by reacting with an aromatic amine and further coupling yield a purple azo dye. Therefore, the color on the test pad changes from light green to purple.
[Specific gravity] Metachromatic. Method
[Leukocytes] Measurement of Leukocyte esterase activity: The test reveals the presence of granulocyte esterases. These esterases cleave an indoxyl ester, and the indoxyl so liberated reacts with a diazonium salt to produce a violet dye.
[pH] pH indicator method: The pH reaction area is impregnated with methyl red and bromothymol blue as indicators. A color change from orange to green indicates pH between 5 and 9.
[Creatinine] Benedict-Behre method²⁾
[Albumin] Protein error of pH indicator: A reaction between tetrabromophenol blue and protein induces a color change from yellow to blue-green.

Parameters: Urobilinogen (URO), Blood (BLD), Red blood cell (RBC), Hemoglobin (Hb), Protein (PRO), Glucose (GLU), Ketones (KET), Bilirubin (BIL), Nitrite (NIT), Specific gravity (S.G), Leukocytes (LEU), pH, Creatinine (CRE), Albumin (ALB).

Components

Reactive ingredients (per test pad 1 cm²)
[Urobilinogen] Metaphosphoric acid: 4.16 mg, 3,4-methylenedioxybenzene diazonium tetrafluoroborate: 0.042 mg
[Blood] Cumen hydroperoxide: 0.223 mg, 3,3',5,5'-tetramethylbenzidine: 0.286 mg
[Protein] Tetrabromophenol blue: 0.015 mg
[Glucose] Glucose oxidase: 0.059 mg, peroxidase: 0.023 mg, 3,3',5,5'-tetramethylbenzidine: 0.088 mg
[Ketones] Glycine: 6.64 mg, sodium nitroprusside: 0.188 mg
[Bilirubin] 2,4-dichlorobenzene diazonium tetrafluoroborate: 0.022 mg
[Nitrite] Sulfanilamide: 0.26 mg, N-(1-Naphthylamino)-3-propanesulfonic acid: 0.076 mg
[Specific gravity] Methylene blue: 0.007 mg, dextran sulfate sodium: 0.016 mg
[Leukocytes] 3-(N-Toluenesulfonyl-L-alanyloxy)-indole: 0.020 mg, 2-methoxy-4-(N-morpholino)-benzenediazonium salt: 0.007 mg
[pH] Methyl red: 0.002 mg, bromothymol blue: 0.038 mg
[Creatinine] 3,5-dinitrobenzoic acid: 0.48 mg
[Albumin] Tetrabromophenol blue: 0.010 mg

Warnings and precautions

- When using this product with Sysmex Urine Chemistry Analyzer UC-1000, please read the analyzer's Instructions for Use carefully.
- Diagnosis should not be made based on the analysis using the MEDITAPE UC-12S test strip alone; use other test results and clinical symptoms for comprehensive diagnosis.
- Use the test strips as directed in this Instructions for Use, otherwise the product quality and safety will not be guaranteed.
- Use disposable powder free gloves to prevent infection during the analysis and disposing the test strips.
- Do not use a product that is suspected to have been frozen.
- Reagents must not be used after its expiry date.
- Avoid refrigeration as much as possible. If the test strip is refrigerated for long-term storage, take it out and wait until it returns to 20 - 25 °C before use.
- Do not use deteriorated, discolored or blackened test strips.
- Do not touch the test pads. Make sure the test strip is free of dirt before use.
- Minimally handle test strips after opening as test strips deteriorate in humid conditions.
- Do not take the desiccant out of the test strip container.
- The test strips are for single use only. Do not reuse.
- Handle samples as biohazardous material.
- Take out a required number of test strips from the container when needed and seal the container immediately after that. Test strips, once taken out, must not be put back into the container. Also, do not replace test strips into other containers.
- Results from visual measurements and analyzer measurements are not always consistent due to some variation between visual recognition and the optical system of the analyzer.

Examination procedure

The tools, apparatus and samples required

A urine collection cup, tissue paper and a method of timekeeping.

The measurement (operation) method

A. Visual measurement

- Following comprehensive mixing of fresh urine, soak test pads completely into fresh urine for a few seconds and then take it out. Dab the test strip with the edge of the cup and tissue paper, etc., to remove excess urine from the test strip.
- Make an assessment of the color change after waiting for the time specified for each parameter in INTERPRETATION OF RESULTS in comparison with the reference color table. Wait a further 60 seconds to make more definite assessment for leukocytes (-) or (I+), in the event the initial color change is not significant enough to make an accurate determination. The waiting time may be 30 to 60 seconds for the screening purpose, however it should be 60 to 120 seconds in an assessment for leukocytes.

B. Analyzer measurement

- Following comprehensive mixing of fresh urine, soak test pads completely into fresh urine for a few seconds and then take it out. Make sure to dab at the test strip with the edge of the cup and tissue paper, etc., to remove excess urine from the test strip.
- Insert the test strip onto the analyzer as specified.
- Upon initial measurement, reflectance is measured on the test strip according to the specified measuring wavelength at the specified time for each measurement parameter to make assessments based on the prearranged respective working curves and output, which then are output as a result.

Storage and shelf life of unopened product

Store the test strips at 1 - 25 °C. Keep away from humidity, direct sunlight and heat. When the product is properly stored in its sealed container, it is stable until the expiration date printed on the label.

Control procedure

Users are responsible for determining the appropriate quality control procedures for their laboratory and for complying with applicable laboratory regulations.
For quality control, use of at least two levels urine control which are commercially available is recommended.

Biological reference intervals

Parameter	Standard range	Parameter	Standard range
URO	0.03-0.97 mg/dL ³⁾	NIT ⁴⁾	
BLD	< 5 cells/HPF ⁴⁾	S.G	1.005-1.030 ¹⁾
PRO	< 30 mg/dL ⁴⁾	LEU	< 12 cells/ μ L ⁴⁾
GLU	2-20 mg/dL ¹⁾	pH	4.5-7.5 ⁴⁾
KET	$\leq$ 2 mg/dL ¹⁾	CRE	0.5-1.5 g/day ¹⁾
BIL	$\leq$ 0.05 mg/dL ⁴⁾	ALB	$\leq$ 23.8 mg/L ¹⁾

Interpretation of results

A. Visual measurement

Parameter		Time	Assessment						
URO		10 sec.	normal	1+ 2.0 (34)	2+ 4.0 (68)	3+ 8.0 (135)	4+ 12.0 (202)	mg/dL (μ mol/L)	
BLD	RBC	30 sec.	-	++ 10	1+ 20	2+ 50	3+ 250	c/ μ L*	
	Hb		-	++ 0.03	1+ 0.06	2+ 0.15	3+ 0.75	mg/dL	
PRO		0 sec.	-	++ 15 (0.15)	1+ 30 (0.3)	2+ 100 (1.0)	3+ 300 (3.0)	4+ 1000 (10)	mg/dL (g/L)
GLU		60 sec.	-	++ 50 (2.8)	1+ 100 (5.6)	2+ 250 (14)	3+ 500 (28)	4+ 2000 (111)	mg/dL (mmol/L)
KET		30 sec.	-		1+ 10 (0.93)	2+ 30 (2.8)	3+ 80 (7.4)		mg/dL (mmol/L)
BIL		20 sec.	-		1+ 0.5 (8.6)	2+ 1.0 (17)	3+ 2.0 (34)		mg/dL (μ mol/L)
NIT		30 sec.	-	0.1	0.3				mg/dL
S.G		30 sec.	1.000	1.005	1.010	1.015	1.020	1.025	1.030
LEU		60 to 120 sec.	-		1+ 25	2+ 75	3+ 500		c/ μ L*
pH		0 sec.		5	6	7	8	9	
CRE		60 sec.		10 (0.1)	50 (0.5)	100 (1.0)	200 (2.0)	300 (3.0)	mg/dL (g/L)
ALB		60 sec.		10 (0.01)	30 (0.03)	80 (0.08)	150 (0.15)		mg/L (g/L)

*c/ μ L = cells/ μ L

Ketones are measured using lithium acetoacetate as standard reference material.

The protein / creatinine ratio (P/C ratio) and the albumin / creatinine ratio (A/C ratio) are automatically computed according to the settings in the analyzer.

P/C ratio (protein / creatinine ratio)

		CRE				
		10 mg/dL	50 mg/dL	100 mg/dL	200 mg/dL	300 mg/dL
PRO	-	dilute * ¹⁾			normal	
	+ -			1+		
	1+					
	2+	2+				
	3+, 4+					

A/C ratio (albumin / creatinine ratio)

		CRE				
		10 mg/dL	50 mg/dL	100 mg/dL	200 mg/dL	300 mg/dL
ALB	10 mg/dL	dilute * ¹⁾			normal	
	30 mg/dL					
	80 mg/dL	2+		1+		
	150 mg/dL					

*1 The "dilute" represents that urine is too dilute for accurate calculation of the P/C ratio and the A/C ratio. Therefore, another urine sample should be obtained for retest.

B. Analyzer measurement

Make assessments in accordance with the prearranged respective working curves and output. The assessment criteria are the same as those in **A. Visual measurement** of the test pads for urobilinogen, blood, protein, glucose, ketones, bilirubin, specific gravity, leukocytes and creatinine. For the test pads for nitrite, pH and albumin, the following assessment method is applied.

Parameter	Assessment							
NIT	-	+						
pH	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5 9.0
ALB	10	30	80	150	over	mg/L		

Nitrite (+), sodium nitrite 0.1 to 0.3 mg/dL.

The P/C ratio and the A/C ratio are automatically computed according to the settings in the analyzer.

Parameter	Assessment						
P/C ratio	dilute	normal	1+ 0.15	1+ 0.30	2+ >=0.50	g/gCr	
A/C ratio* ²⁾	dilute	normal	1+ 30	1+ 80	1+ 150	>=1+ >=80	>=1+ >=150 >=300 mg/gCr

*2 For A/C ratio, when albumin is "over" and creatinine is 200 or 300 mg/dL, the qualitative value and the semiquantitative value shall be computed as >=1+ and >=150 or >=80, respectively.

[Note 1] Urobilinogen (-) cannot be confirmed with this method⁴⁾.

[Note 2] Bacteriuria cannot be denied even if nitrite is negative because some bacteria do not reduce nitrate, and in a case where urine lacks nitrate, the test pad will show nitrite negative since nitrate-reducing bacteria cannot produce nitrite⁴⁾.

[Note 3] The test pad for leukocytes measures esterase activity in leukocytes. Therefore, an assessment may be different from a result from urinary sediment depending on extent of disintegration of leukocytes in urine⁵⁾.

[Note 4] The P/C ratios and A/C ratios obtained with the test strips are used for screening purpose, therefore a quantitative measurement based on appropriate treatment guidelines should be performed for diagnosis.

Performance characteristics

1. Sensitivity

- For the test pads for urobilinogen, blood, protein, glucose, ketones, bilirubin, nitrite, leukocytes, creatinine and albumin, when the following 2 reference urine levels are measured, the resulted assessment is consistent with the preset assessment grades and then may be clearly differentiated.
- For the test pad for specific gravity, when a test is performed by using reference urine at 1.000, 1.005, 1.010, 1.015, 1.020, 1.025 and 1.030, the resulted values as outcome are those in the preset assessment grade $\pm$ 0.005.
- The test pad for pH shows outcome in the following manner.

A. Visual measurement

The test pads, when tests are performed by using the reference urine with pH 5, 6, 7, 8 or 9, show values consistent with the preset reference color table as outcome.

B. Analyzer measurement

For the test pad for pH, when a test is performed by using reference urine with pH at 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5 or 9.0, the resulting values as outcome are those in the preset assessment grade $\pm$ 0.5.

Parameter		Reference urine level	Parameter	Reference urine level
URO		≤ 0.2 mg/dL and 2.0 mg/dL	BIL	0 mg/dL and 0.5 mg/dL
RBC	RBC	0 cells/μL and 10 cells/μL	NIT	0 mg/dL and 0.1 mg/dL
	Hb	0 mg/dL and 0.03 mg/dL	LEU	0 cells/μL and 25 cells/μL
PRO		0 mg/dL and 15 mg/dL	CRE	10 mg/dL and 50 mg/dL
GLU		≤ 10 mg/dL and 50 mg/dL	ALB	10 mg/L and 30 mg/L
KET		0 mg/dL and 10 mg/dL		

Ketones are measured using lithium acetoacetate as standard reference material.

2. Accuracy

- For the test pads for urobilinogen, blood, protein, glucose, ketones, bilirubin, nitrite, leukocytes and albumin, when a test is performed by using reference urine in which levels are equivalent to those in each assessment grade, outcome is consistent with the preset assessment grades.
- The test pads for specific gravity and pH show outcome in the same manner as shown in **1. Sensitivity**.
- The test pad for creatinine shows outcome in the following manner.

A. Visual measurement

The test pads, when tests are performed by using reference urine in which levels are equivalent to those in each assessment grade, show outcomes consistent with the preset reference color table.

B. Analyzer measurement

The test pads, when tests are performed by using reference urine in which levels are equivalent to those in each assessment grade, show outcomes consistent with the preset assessment grades in cases of 10, 50 or 100 mg/dL. The assessment grade may become higher by 1 grade than the preset one in cases of 200 mg/dL, and the assessment grade may become lower by 1 grade than the preset one in cases of 300 mg/dL.

3. Repeatability

The test pads show the following outcome, when reference urine with known levels is tested simultaneously for 5 times.

- The test pads for urobilinogen, blood, protein, glucose, ketones, bilirubin, nitrite, leukocytes and albumin show a consistent outcome.
- The test pad for specific gravity shows values in the preset assessment grade $\pm$ 0.005 as outcome.
- The test pad for pH shows a consistent outcome in visual measurement and shows values in the preset assessment grade $\pm$ 0.5 as outcome in analyzer measurement.
- The test pad for creatinine shows a consistent outcome in visual measurement, and, in analyzer measurement, shows a consistent outcome at lower levels (10 to 50 mg/dL) and shows 200 mg/dL or 300 mg/dL at higher levels (200 to 300 mg/dL).

4. Measuring range

Parameter		Measuring range	Parameter	Measuring range
URO		2.0-12.0 mg/dL	NIT	0.1-0.3 mg/dL
BLD	RBC	10-250 cells/ μ L	S.G	1.000-1.030
	Hb	0.03-0.75 mg/dL	LEU	25-500 cells/ μ L
PRO		15-1000 mg/dL	pH	5.0-9.0
GLU		50-2000 mg/dL	CRE	10-300 mg/dL
KET		10-80 mg/dL	ALB	10-150 mg/L
BIL		0.5-2.0 mg/dL		

Ketones are measured using lithium acetoacetate as standard reference material.

5. Correlation

Parameter	Number of patients	Concordance rate (%)
URO	254	99.2
BLD	279	97.1
PRO	224	96.4
GLU	286	97.9
KET	281	99.3
BIL	299	99.0
NIT	288	99.3
S.G	74	97.3
LEU	265	96.2
pH	215	91.6
CRE	88	94.3
ALB	108	90.7

6. Interferences^{4), 6)}

Parameter	Interfering substance
URO	• Not affected by porphobilinogen, indole, para-aminosalicylic acid, sulfonamide ⁹⁾ and urea ⁷⁾ that react to Ehrlich's aldehyde reagent. • The test may show a greenish color in the presence of strongly positive bilirubin urine.
BLD ⁹⁾	• The test may be false negative in the presence of a large amount of reducing agent such as ascorbic acid and nitrite in the urine. The in-house study showed that sodium nitrite concentrations $\leq$ 10 mg/dL displayed no false negative results when the measured hemoglobin concentration was 0.06 mg/dL and 0.15 mg/dL. • The test may be false positive if affected by oxidizing agents such as hypochlorous acid and bleaching powder. The in-house test was false positive in the presence of sodium hypochlorite $\geq$ 1.2 mg/dL. • The reactivity may be inhibited in baruria. • Reaction to myoglobin may occur ⁹⁾ . • The test may be false positive if drugs containing -SH groups (e.g., glutathione agents and buccillamine, etc.) have been taken.
PRO ¹⁰⁾	• The test may be false positive in urine of pH $\geq$ 8 or urine having strong buffering action. • The test may be false positive if detergent or disinfectant (quaternary ammonium compound or chlorhexidine) is left in the container. • The reaction to globulin and mucoprotein, etc. is weaker compared with albumin. • It has been reported that although the reactivity of Bence Jones protein is lower than that of albumin, they often exhibit similar sensitivity ¹⁰⁾ . • The test may be positive in urine containing semen.
GLU ¹²⁾	• False negative test results are possible in the presence of high levels of ascorbic acid. Investigations found that ascorbic acid concentrations $\leq$ 200 mg/dL display no false negative results when the measured glucose concentration is 100 mg/dL. • The in-house test showed no effect of sodium chloride (3 %), uric acid (150 mg/dL) or sodium nitrite (10 mg/dL). • The test may be false positive if affected by oxidizing agents such as hypochlorous acid and bleaching powder. The in-house test was false positive in the presence of sodium hypochlorite $\geq$ 1.2 mg/dL. • The test pad reacts to galactose. • The reactivity may be inhibited in baruria.
KET	• The test may be false positive or an abnormal color may be shown in the presence of a large amount of phenylpyruvic acid, pyruvic acid, oxaloacetate, α -ketoglutaric acid or phenolsulfonphthalein (PSP) ¹⁰⁾ . • The test pad does not react to β -hydroxybutyric acid. • The test pad may be false positive if drugs containing -SH groups (e.g., glutathione agents and buccillamine) have been taken ¹⁰⁾ .
BIL	• The test may be false negative in the presence of a large amount of ascorbic acid or nitrite. • The test may be false positive in the presence of a large amount of urobilinogen or 5-hydroxyindoleacetic acid (5-HIAA) ¹⁰⁾ . • When etodolac has been taken, the test may be false positive showing pink (unlike the usual color of bilirubin) due to the reaction to its metabolite phenol derivative ¹⁰⁾ .
NIT	• The test may be false negative in the presence of a large amount of ascorbic acid. The in-house study showed that ascorbic acid concentrations $\leq$ 25 mg/dL displayed no false negative results when the measured nitrite concentration was 0.1 mg/dL.
S.G ¹⁶⁾	• Not react with non-ionic substances such as glucose in urine and urea, and with proteins. • Not affected by the pH of urine.
LEU ¹⁷⁾	• The test may be false positive in the presence of formaldehyde (urine preservative) ⁹⁾ . • The test may be false negative under the presence of protein $\geq$ 500 mg/dL. • The test may be false negative under the presence of cephalixin, gentamicin or boric acid (urine preservative). • The test may show a different color reaction than the reference color in the urine containing bilirubin or nitrofurantoin.
pH	• Use of volatile materials such as acid and alkali may affect the test.
CRE	• The assessment value may be lower in urine containing ketones. However, the in-house test showed no effect of lithium acetoacetate $\leq$ 20 mg/dL.
ALB ¹⁸⁾	• The assessment value may be higher in urine of pH $\geq$ 8, urine having strong buffering action, or urine containing alkaline mucus. • The test may be false positive if detergent or disinfectant (a quaternary ammonium compound or chlorhexidine) remains in the container. • The test may be false positive in the presence of a large amount of hemoglobin or myoglobin or macroscopically-obvious hematuria. The test may be positive in urine containing semen. • With regard to the following proteins, no impact was observed up to the concentration level shown in parentheses ¹⁾ : Lysozyme (50 mg/L), α 1-microglobulin (25 mg/L), β 2-microglobulin (15 mg/L), prealbumin (200 mg/L), immunoglobulin (500 mg/L), α 1-acidic glycoprotein (300 mg/L), haptoglobin (25 mg/L), transferrin (25 mg/L), α 1-antitrypsin (25 mg/L), Bence Jones protein (100 mg/L) and retinol-binding protein (15 mg/L)

Limitations of the examination procedure

Severely discolored urine and urine discolored by medication may affect the analysis because of abnormally colored test pads.^{6), 19)}

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